
1. Other functionalities


Welcome window:

To enter the system, you should click on the button with the arrow.




Login:

To log in the system, you have to be registered and enter your email and password. In case you are not allowed it could be because your institution is suspended or not approved yet.




Register an institution and its manager:

If you are not registered and you want to be, you can sign up in this page. You can access it from the login page in the "Do you want to register an institution?" link. Until your institution is not approved, you cannot login in the application.



Home page:

This is the home page. Here you can also see the sidebar, that can be folded and unfolded. From the sidebar you can access the other parts of the application.




View the profile:

By clicking the first button of the lower group in the sidebar you can see your profile and change your password.




Change language:

In the indicated button you can change the language from English to Spanish.



About and credits:

If you click the about button you can see a brief description of the application and project and the credits document.



Help:

You can access this manual by clicking on the help button.




Logout:

You can log out by clicking on this button in the sidebar. It will redirect you to the login page.



Search:

You can filter all the lists using the search bar.




2. Users

You can only access user-related content if you are a manager or an administrator.


Users:

In the users list you can click a user to display its information.




Create a user:

By clicking the button at the right of the search bar you can add a user.




Modify a user:

At the right of a user's name, you can open a drop down, and modify users.



Delete a user:

At the right of a user's name, you can open a drop down, and delete users.




3. Institutions

You can only access institution-related content if you are an administrator.


Institutions:

In the institution list, accessible from the sidebar, you can click on an institution to display its information.




Create an institution:

In the "+" button you can create institutions and their managers.




Modify an institution:

You can modify an institution by clicking on the button with three dots at the right of the institution's name.




Suspend an institution:

You can suspend an institution by clicking on the button with three dots at the right of the institution's name. Its users won't be able to access the application from then on.



Pending institutions:

If you select the pending review filter, institutions with a pending status will appear. If you click on the institution, you can approve or dismiss their request to be in the system. Until you approve an institution, their users will not be able to access the application.




4. Subjects

You can only access subject-related content if you are a technician or manager.


Subjects:

In the institution list, accessible from the sidebar, you can click on an institution to display its information.




Subject's page and evaluate a subject:

When clicking on a subject of the list, the system redirects you to its personal page, where you can check its information and click on the evaluate button. For the system to evaluate a subject correctly, you should have uploaded first the files.




Upload tests:

You can upload the eye tracker data to the application in the button at the right of "Tests". You should upload the interchanges file and the fixations and time files. The files must be in a .tsv format.



Move a subject:

You can move a subject by clicking on the three dots at the right of "Personal information". There you can select to which institution you want to move the subject.



Confirm a diagnosis:

After having evaluated a subject, you can confirm their diagnosis in the "Add diagnosis confirmation" button, which is not accessible until the automatic result is saved.

